



MAKERS OF RESPIRATORY VALVES SINCE 1938


HANS RUDOLPH, inc.

www.rudolphkc.com

TEL: (913) 422-7788 U.S.A. & CANADA (800) 456-6695
 FAX: (913) 422-3337 E-Mail: hri@rudolphkc.com
 8325 COLE PARKWAY, SHAWNEE, KANSAS 66227 U.S.A.

Data Sheet...

MOUTHPIECES, NOSE CLIPS, SAMPLING PORTS, COUPLERS, TUBING & ACCESSORIES



Mouthpiece Saliva Trap Type (Green)			
<i>Reusable Series 9060</i>			
Silicone Rubber Maintenance Document: 691207			
Part Number	Bite Size	Port Size Range Diameter (mm)	Dead Space (ml)
602077	Large	35-38	93
602078	Medium	28-34	66
602079	Small	22-27	44




Mouthpiece Standard Type (Clear)			
<i>Reusable Series 9060</i>			
Silicone Rubber Maintenance Document: 691207			
Part Number	Bite Size	Port Size Range Diameter (mm)	Dead Space (ml)
602070	Large	35-38	39.5
602073	Medium	28-34	23
602076	Small	22-27	12.5

Nose Clip Reusable Series 9015	
Styrenic Thermoplastic Elastomer, Neoprene Foam, Stainless Steel	
Maintenance Document: 691207	
Part Number	Unit of Measure
112571	Each

Nose Clip Disposable Series 9014	
Polypropylene Plastic, EVA Foam	
Part Number	Unit of Measure
201413	50 piece bag
201414	2000 piece bag

Sampling Ports (Hose Barb & Female Luer Types)				
Nylon Plastic				
Part Number	Hose Barb (B) Female Luer (L)	Barb Fits Tubing ID mm(in)	Plug or Cap P/N	Installation Part Number
171161	B	1.6 (.062)	669023	200120
171095	B	2.4 (.093)	669023	200121
171163	B	3.2 (.125)	669021	200122
171179	L	Male Luer	171180	200123



Couplers (Male or Female Luer with opposite end Hose Barb)		
Nylon Plastic		
Part Number	Male Luer (M) Female Luer (F)	Barb Fits Tubing ID mm(in)
171210	M	1.6 (.062)
171211	M	2.4 (.093)
171212	M	3.2 (.125)
171269	F	1.6 (.062)
171270	F	2.4 (.093)
171271	F	3.2 (.125)



G



Tubing Flexible (Pneumatics and gas sampling lines) **G**
Tygon, Polyurethane, PVC

Part Number	Material	Size OD/ID mm (in)	Color	Pressure Rating Bar (PSI)
666060	Tygon	3.2(.125) / 1.6(.062)	Clear	2.75 (40)
666061	Tygon	4.0(.156) / 3.2(.125)	Clear	2.75 (40)
666062	Tygon	4.8(.187) / 3.2(.125)	Clear	2.75 (40)
666064	Tygon	7.9(.312) / 8.4(.187)	Clear	2.75 (40)
666065	Tygon	9.5(.375) / 6.4(.250)	Clear	2.75 (40)
666014	Polyurethane	3.2(.125) / 1.6(.062)	Clear	6.9 (100)
666025	Polyurethane	6.4(.250) / 3.2(.125)	Blue	6.9 (100)
666016	Polyurethane	9.5(.375) / 6.4(.250)	Gray	6.9 (100)
666015	PVC	6.4(.250) / 3.2(.125)	Clear	4.69 (68)

Accessories Miscellaneous **H**

Item # See Photo	Part Number	Description
1	990001	Calipers 6"
2	990002	Allen Wrench Set 12 Piece
3	901000	Drill & Tap Set #6-32
3	901001	Drill & Tap Set #8-32
3	901002	Drill & Tap Set #10-32
4	171218	Retaining Ring Pliers (Ballon Valves)
5	905000	Oring Removal Tool (pic)
6	905001	Wrench 1/4" Box & Open End
7	667009	Adhesive Cyanoacrylate 1 oz
8	660170	Silicone Lubricant 5.3 oz
9	871002	Teflon Tape Roll 1/4" width (threads)
10	602032	Component Grip Pad



H